



**PAIN SCORES**

0 = No pain

The patient Total = Volume infused + bolus given.

Please check the pump strictly every 8 hours for bolus attempts, bolus given, and patient total. To do this:

1. Press the red stop/no button.
2. Press the blue info key to obtain the **8 HOURLY SCREEN. THIS IS NOT A RUNNING TOTAL.**
3. Press the stop/no button twice, then press the green start/ok button to resume the infusion.
4. Check VTBI (remaining volume) by pressing the info key twice.
5. At the same time each day, complete a 24-hour patient total by following step 4, and complete the following page.

1= Mild Discomfort

2= Moderate Discomfort

3= Painful

4=Severe Pain

5= Excruciating

**24 HOUR PATIENT TOTAL SUMMARY** (At the same time each day)  
This is a cumulative total, therefore you will need to subtract the previous 24-hour total

Date	Time	Patient Total	Date	Time	Patient Total

Date and Signature									
<b>Medication Bag Change</b>									
<b>Line / Filter Change</b>									
<b>Dressing changed</b>									

**Dressing Changes weekly- See Action Card 18 or 'Intrathecal Analgesia in Cancer Pain: A resource for Health Professionals'.**